



The quality
your patients' lungs
deserve.

The new PARI BOY® Go –
the quietest and most mobile
PARI BOY ever.


PARI BOY® Go

Premium quality – easy to use

Inhalation therapy depends on consistent use. Yet when treatment is perceived as cumbersome or time-consuming, adherence may suffer.

The new PARI BOY Go is a high-quality, mobile addition to the PARI BOY product line. It supports your patients in their inhalation therapy with its intuitive use and ergonomic design.

Thanks to its compact size, the PARI BOY Go can be used flexibly in everyday life and on the go, while mask use additionally allows inhalation in a lying position.



PARI care mode

for reliable performance

With its membrane nebuliser technology, the PARI BOY Go is particularly quiet and powerful. The membrane is the heart of the device.

The PARI BOY Go reliably reminds your patient to maintain the membrane at regular intervals with the convenient PARI care mode – completely automatically. It ensures consistently high performance throughout the nebuliser's service life by backflushing the membrane to remove any residual inhalation solution from the membrane pores.



For a video on easy handling,
cleaning and care, visit
www.youtube.com/parigmbh



PARI BOY® Go

at a glance

Which indication?

Acute and chronic respiratory conditions such as bronchitis, asthma, and COPD

Who is it for?

Approved for children 2 years and older and adults

For multi-patient use?

Controller: multi-patient use
Nebuliser top: single-patient use

Used with?

For all medications approved for inhalation, except for inhaled antibiotics



Excellent aerosol performance

MMAD: ~4.4 μ m
RDDR: ~125 μ l/min
RF: ~58 %
Residual volume: ~0.4 ml

Package Contents



Controller
(incl. batteries)



Nebuliser
top



Mouthpiece



Child mask*



Adult mask*



Bag

The nebuliser top and mouthpiece are available as a Year Pack for convenient replacement.

Compelling benefits

for you and your patients



Small and mobile
wherever you Go!



Virtually silent
for discreet nebulisation



Can be used while lying down
and using a mask



Dual power supply
(battery + USB-C)
ready to use at any time



Short inhalation time
approx. 4 min. for 2.5 ml NaCl



Thermal disinfection
(boiling)

* Please note that the included package contents may vary by country.

You can count on PARI

PARI – a trusted brand.

With over 100 years of experience, we develop and manufacture our inhalation devices at our facilities in Upper Bavaria, Germany, and we are committed to the highest quality – both in our products and in our customer service.

As a trusted partner in the healthcare industry, we set the standard:

- Made in Germany – quality, durability, and reliability
- 1968: PARI BOY® set the milestone in home nebuliser therapy
- In-house development, manufacturing, and quality assurance
- Long-term availability of spare parts
- Reliable after-sales service



We're here to help!

If you need information about a repair or service for a PARI device, please visit www.pari.com/partners to find a distributor and/or service partner in your country.

We look forward to hearing from you. Simply give us a call or send an e-mail to export@pari.de for a prompt response.

PARI GmbH
Moosstr. 3, 82319 Starnberg, GERMANY
phone: +49 (0) 8151 279 220
fax: +49 (0) 8151 279 6220
e-mail: export@pari.de
www.pari.com

Distributed by: