



PARI CENTRAL AIR

Technical data sheet

Pressure and flow reducer for use with a PARI nebuliser in professional health institutions.

PARI CENTRAL AIR reduces the gas pressure and flow of the central medical gas supply system for an optimum function with a connected PARI nebuliser. It is used instead of a compressor. Connectors can be used for the central gas supply according to following standards: DIN 13260-2 and SS 87 524 30



Example for DIN 13260-2 (connector varies between standards)

Contents:

- PARI CENTRAL AIR
- Connection tubing long (f/f) with tubing adapter (Length of tubing: 1.9m)
- PARI Central Clip (not included in version SS 87 524 30)

Technical Characteristics:

	Minimum	Nominal	Maximum
Inlet pressure:	2.6 bar/260 kPa	5.0 bar/500 kPa	6.0 bar/600 kPa
Outlet/operating pressure:	0.6 bar/60 kPa	1.6 bar/160 kPa	1.9 bar/190 kPa
Flow:	3.0 l/min	5.0 l/min	6.0 l/min

Dimensions:

	Length	Diameter	Weight
PARI CENTRAL AIR (DIN 13260-2)	86 mm	26 mm	67 g
PARI CENTRAL AIR (SS 87 524 30)	99 mm	26 mm	80 g

PARI CENTRAL AIR

- does not require any electricity to operate
- needs no time-consuming adjustments or calibration, just plug it in
- its operating flow and pressure are pre-set
- for use with PARI jet nebulisers (except PARI LC SPRINT SINUS) and all suitable accessories
- is made from high-quality materials (aluminium and stainless steel) and is maintenance-free.
- can be easily disinfected in accordance with instruction in IFU
- comes with the PARI CENTRAL Clip – mount for the nebuliser and / or tubing whilst filling the nebuliser or when not in use

Overview of Item Numbers

Item Number	Item Description	IFU Languages
058G1100	PARI CENTRAL AIR – DIN 13260-2	de, en, fr, it, nl
058G1110	PARI CENTRAL AIR – SS 8752430	da, fi, no, sv

Overview of Nebuliser Hospital Packs

Item Number	Item Description	IFU Languages
023G5800	PARI LC SPRINT CENTRAL	da, de, en, nl, no, sv
023G6010	PARI LC SPRINT Hospital Pack (pack of 10)	da, de, en, nl, no, sv
023G6011	PARI LC SPRINT Junior Hospital Pack (pack of 5)	da, de, en, nl, no, sv
023G6012	PARI LC SPRINT STAR Hospital Pack (pack of 5)	da, de, en, nl, no, sv

All PARI Jet Nebulisers (excl. PARI LC SPRINT SINUS) are suitable for use with PARI CENTRAL

Spare Parts

Item Number	Item Description	IFU Languages
058G2001	PARI CENTRAL Clip (Pack of 5)	n.a.
041G4550	PARI Connection Tubing for PARI CENTRAL (with tubing adapter)	n.a.



PARI GmbH
Moosstraße 3, 82319 Starnberg, GERMANY
E-Mail: export@pari.de
Phone: +49 (0) 8151 279 220
www.pari.com